# QUICK REFERENCE GUIDE

Easy operation — fast results

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---

 Direct current.

 Caution. Refer to any accompanying documents.

 Biological risks. Indicates that there are potential biological risks associated with the medical device. See “Universal Precautions” on pages 1-2.

 USB connection

 Network connection

 Manufacturer

 Date of manufacturing in the USA

 Authorized Representative of the European Union

 Importer contact (EU)

 In vitro diagnostic medical device

 Serial Number

 Catalog Number

 Device may be used for Near-Patient-Testing

 Do not dispose as household waste. Separate collection and recycling.

 Unique Device Identifier in human and machine-readable form used to adequately identify medical devices through their distribution and use.

 UDI Barcode structure in Health Industry Bar Code (HIBC) standard

 Use by

 Batch code

 Consult instructions for use

 X number of Test Devices in kit

 Manufacturing Sequence

 Temperature Limitation

 Part Number

 Denotes conformity to specified European directives and regulations

This product is designed as an In Vitro Diagnostic (IVD) system in accordance to the European Directive on In Vitro Diagnostic Medical Devices (2017/746/EU).

## Index

---

|                             |    |
|-----------------------------|----|
| Quick Reference Guide       | 5  |
| Kurzanleitung               | 19 |
| Guide de référence rapide   | 35 |
| Guida rapida di riferimento | 51 |
| Guía rápida de referencia   | 67 |
| Guia rápido de referências  | 83 |
| Ghid de referință rapid     | 99 |

English

Deutsch

Français

Italiano

Español

Português

Română

## Intended Use

The *Piccolo Xpress*® chemistry analyzer provides quantitative in-vitro determinations of clinical chemistry analytes in lithium-heparinized whole blood, heparinized plasma, or serum. The *Piccolo Xpress*® will be referred to as the “analyzer” throughout this guide.

**CAUTION:** The analyzer may not operate as intended if used other than described in this guide. It may produce inaccurate or no results, and may pose a safety hazard.

**Note:** *Use only Piccolo reagent discs with this analyzer.*

### Applicable to U.S. Customers Only:


CLIA WAIVED: The Lipid Panel, Lipid Panel Plus, Liver Panel Plus, General Chemistry 6, General Chemistry 13, Kidney Check, Comprehensive Metabolic Panel, Basic Metabolic Panel, Electrolyte Panel, MetLyte 8, and Renal Function Panel Reagent Discs are CLIA Waived. These discs are waived under CLIA '88 regulations. For CLIA waived labs, only lithium heparin whole blood may be tested. For use in moderate complexity labs, lithium heparinized whole blood, lithium heparinized plasma, or serum may be used. A CLIA Certificate of Waiver is needed to perform CLIA waived testing. A Certificate of Waiver can be obtained from the Centers for Medicare & Medicaid Services (CMS). If a waived laboratory modifies the test system instructions, then the test is considered high complexity and subject to all CLIA requirements.

### Universal Precautions

Operator health and safety regulations require that Universal Precautions be observed at all times while handling human blood samples or working with this analyzer.

# Unpacking and Setup

- 1 Remove the analyzer from the shipping carton.
- 2 Place it on a level surface that is free of hair, dust, and other contaminants. Do not locate it near a sunny window or any other heat source. Make sure the analyzer is at least 20 cm away from any wall to provide adequate ventilation and access to the power connection, Ethernet, and USB ports.
- 3 Check the components received with the *Piccolo Xpress*® per the figure to make sure everything required to set up the analyzer is included.

**Note:** *The shipping carton also includes following items: USB cable, warranty card, and accessory kit.*

- 4 Complete the warranty card and mail it to Abaxis within 10 days to start the warranty period.
- 5 Download the Piccolo Xpress Operator's Manual here:
- 6 Open the printer cover.
- 7 Unwrap the paper roll and unroll several centimeters.
- 8 Put the paper into the printer so that it unrolls from the bottom toward the front of the analyzer as shown. Make sure several centimeters of paper extend out of the printer slot.
- 9 Press the cover closed until it locks into place.
- 10 Pull gently on the end of the paper until it is taut.

- 11 Plug the power jack into the analyzer, then plug the detachable power cord into the power adapter and into a grounded electrical outlet.

***Note:** To prevent power surges or drain, do not plug the analyzer into the same circuit as a centrifuge or any other high-current device. Abaxis also recommends using a suitable surge protector.*

- 12 Press the power button to turn on the analyzer. The display will show:


- 13 After passing the self-test and reaching operating temperature, the analyzer is ready to run the first reagent disc. The display will read **Analyze**. The analyzer may require up to 4 minutes to reach operating temperature.
- 14 Check the analyzer date and time to be sure they are correct.
- 15 Reference ranges are preset by the factory.

**Please make sure to be familiar with the system before running any patient samples.**

## Filling the Reagent Disc

A 100  $\mu\text{L}$  micropipette is provided with the analyzer.


- 1 Take a reagent disc out of its pouch. Discs can be used directly from the refrigerator without warming.
- 2 Attach a new tip to the end of the pipette.
- 3 Fill the pipette with the sample. Make sure there are no air bubbles or gaps.
- 4 Tilt the disc to  $45^\circ$  with the sample port above the fill line, so that the entire sample will flow into the sample chamber.
- 5 Place the pipette tip into the disc sample chamber and release the sample. Make sure to keep the pipette plunger pressed down until the tip is removed from the sample port.
- 6 Discard the pipette tip in a biohazard container.
- 7 Clean any sample spilled on the disc with a lint-free tissue. Take care that the tissue does not withdraw any sample from the sample port. Dispose of the tissue in a biohazard container.

In addition to the micropipette provided, any other transfer device can be used to add approximately 100  $\mu\text{L}$  of sample into the sample chamber.

*Note: When using lithium heparin whole blood tubes, gently rock them several times before use.*

# Running a Patient Sample

- 1 Turn on the analyzer by pressing the Power button on the front. The analyzer will start up and perform a self-test. **Warming** is displayed during the warm-up time. **Analyze** appears when the analyzer reaches operating temperature.
- 2 Select **Analyze** to open the disc drawer.
- 3 Place the disc in the drawer.

- 4 Select **Close** to close the disc drawer.
- 5 Select **Patient** or **Control** to run a patient sample or serum control.
- 6 Enter an ID number for the sample (up to 14 characters), then select **Done**.
- 7 The analyzer checks the disc type and begins processing the sample.
- 8 When the sample is processed, the analyzer stores the results and shows that the analysis is complete.
- 9 Select **Open** to open the drawer.
- 10 Discard the used disc in a biohazard container.
- 11 To analyze another sample, insert a new disc and repeat the above procedure.
- 12 When finished, select **Close** to return to the home screen.

**Note:** The Piccolo Xpress® chemistry analyzer accepts lithium-heparinized whole blood, plasma, or serum for patient samples. Lithium heparin is the only anticoagulant recommended for use with the Piccolo Xpress®.

# Interpreting Results

## Test results:

- Analyte
- Analyte Concentration
- Reference range
- Units

## QC OK:

iQC® checks passed

## Head of result printout:

- Disc type
- Date/Time
- Sample Type
- Sample ID
- Lot #
- Serial #

| Piccolo Xpress     |           |           |            |  |
|--------------------|-----------|-----------|------------|--|
| Comprehensive      | Metabolic |           | Panel      |  |
| 11 Apr 2018        |           |           | 01:15 PM   |  |
| Sample Type:       |           |           | Patient    |  |
| Patient ID:        |           |           | 77752      |  |
| Disc Lot Number:   |           |           | 8061BB4    |  |
| Serial Number:     |           |           | 0000P20587 |  |
| .....              |           |           |            |  |
| NA+                | 137       | 128-145   | mmol/L     |  |
| K+                 | 4.0       | 3.6-5.1   | mmol/L     |  |
| tCO2               | 28        | 18-33     | mmol/L     |  |
| CL-                | 103       | 98-108    | mmol/L     |  |
| GLU                | 5.3       | 4.05-6.55 | mmol/L     |  |
| CA                 | 2.2       | 2.0-2.58  | mmol/L     |  |
| BUN                | 4.2       | 2.5-7.9   | mmol/L     |  |
| CRE                | 79        | 53-106    | umol/L     |  |
| eGFR               | 60        | c         | mL/min     |  |
| ALP                | 91        | 42-141    | U/L        |  |
| ALT                | 29        | 10-47     | U/L        |  |
| AST                | 25        | 11-38     | U/L        |  |
| TBIL               | 14.2      | 3.4-27.4  | umol/L     |  |
| ALB                | 44        | 33-55     | g/L        |  |
| TP                 | 78        | 64-81     | g/L        |  |
| .....              |           |           |            |  |
| QC OK              |           |           |            |  |
| HEM 1+ LIP 0 ICT 0 |           |           |            |  |

## HEM/LIP/ICT:

Qualitative indices  
for Sample interferences

## Symbols printed on the printout

- \* result outside the reference range
- < result below the dynamic range of the analyzer
- > result above the dynamic range of the analyzer
- c calculated value

**QC OK** all internal quality control (iQC®) checks were passed (see below)

Qualitative indices show the degree of hemolysis (HEM), icterus (ICT), and lipemia (LIP) found in the sample on a scale of: 0 (clear), 1+ (slight), 2+ (moderate), and 3+ (gross).

The analyzer features an internal quality control (iQC®) system that suppresses values and displays warning or error messages when a problem occurs. The iQC® conducts several quality and sufficiency checks with every run, while monitoring all technical components of the analyzer. A suppression or warning does not always mean that the analyzer has a technical fault. It proves that the iQC® is working properly and ensures that only accurate results will be released.

Refer to “Troubleshooting” Section for a list of error messages.

- Warnings and error messages include an internal error code that will assist the Abaxis Technical Support in diagnosing the issue. Record the error message and/or print the error report before contacting Abaxis Technical Support.

#### *Customer and Technical Support:*

USA/Canada: 1-800-822-2947 or [ab-piccolotechsupport@zoetis.com](mailto:ab-piccolotechsupport@zoetis.com)

International: +49 (6155) 780 210 or [int.piccolosupport@zoetis.com](mailto:int.piccolosupport@zoetis.com)

# Quality Control

## International Quality Control:

Analyzer performance can be verified by running controls. Please contact Abaxis Technical Support for a full list of approved quality control materials with acceptance ranges.

## The following is applicable to US Quality Control materials only.

Bioresource Technology assay controls (BRT Controls) are recommended. Other human serum or plasma-based controls may not be compatible. Quality control materials should be stored per their package insert.

BRT controls and number of control levels recommended:

| Type of Control                           | Number of Control Levels | Storage Condition                                                                | Thawed and Unopened                                                                                                       | Opened                                                                                                       |
|-------------------------------------------|--------------------------|----------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------|
| Verification Samples Liquid               | 3 levels                 | <b>Unopened</b> at $\leq -15^{\circ}\text{C}$ . Store vials away from the light. | Control can be used for up to 14 days when stored <b>unopened</b> at $28^{\circ}\text{C}$ . For optimum Bilirubin and CO. | Cap and store at $2-8^{\circ}\text{C}$ . Dispose of at the end of day or upon completion of data collection. |
| Chemistry Control Liquid Assayed          | 2 levels                 | <b>Unopened</b> at $\leq -15^{\circ}\text{C}$ . Store vials away from the light. | Control can be used for up to 14 days when stored <b>unopened</b> at $28^{\circ}\text{C}$ . For optimum Bilirubin and CO. | Cap and store at $2-8^{\circ}\text{C}$ . Dispose of at the end of day or upon completion of data collection. |
| Verification Samples with CRP Liquid      | 3 levels                 | <b>Unopened</b> at $\leq -15^{\circ}\text{C}$ . Store vials away from the light. | Control can be used for up to 14 days when stored <b>unopened</b> at $28^{\circ}\text{C}$ . For optimum Bilirubin and CO. | Cap and store at $2-8^{\circ}\text{C}$ . Dispose of at the end of day or upon completion of data collection. |
| Chemistry Control with CRP Liquid Assayed | 2 levels                 | <b>Unopened</b> at $\leq -15^{\circ}\text{C}$ . Store vials away from the light. | Control can be used for up to 14 days when stored <b>unopened</b> at $28^{\circ}\text{C}$ . For optimum Bilirubin and CO. | Cap and store at $2-8^{\circ}\text{C}$ . Dispose of at the end of day or upon completion of data collection. |

## Procedure

The control samples should be treated the same as a patient sample and run according to the instructions with the instrument, kit, or reagent being used. Before sampling, the control samples should be mixed thoroughly but gently.

Thaw BRT samples at room temperature (18–25°C) for 1 hour or until completely thawed. Mix the vial by rocking it several times. Before sampling, gently swirl until mixed completely, with no visible signs of solid material. Do not shake the vial.

After sampling, promptly recap the control samples and store at 2–8°C. Dispose of the control samples at the end of day or when collection is finished.

## Results

The mean values and acceptable ranges are lot-specific for each control. Individual laboratory values should fall within the acceptable ranges identified in the control leaflet.

If control results are out of range, repeat one time. If still out of range, call Abaxis Technical Support. Do not report results if controls are outside their labeled limits. See the *Operator's Manual* for a detailed discussion on running, recording, interpreting, and plotting control results.


Abaxis recommends control testing as follows:

- at least every 30 days
- whenever laboratory conditions have changed significantly
- when training or retraining of personnel is indicated
- when test results do not match patient symptoms or clinical findings
- with each new lot (CLIA waived tests in waived status labs)

## Cleaning and Maintenance

The analyzer requires minimal maintenance. Regular care of the exterior, screen, and air filter will help ensure reliable operation – frequency will depend on environmental conditions.

- Clean the exterior and case weekly with mild detergent and a soft, damp cloth.
- Clean the screen periodically using a soft, lint-free cloth dampened with a glass-cleaning fluid or window cleaner. Disinfect the screen using a 10% bleach solution (1 part bleach to 9 parts water) applied to a lint-free cloth. Do not use any cleaner containing alcohol. Do not spray cleaner directly onto the display – dampen the cloth instead.
- Clean air filter twice per year. Check it more often if the analyzer is in an environment with a lot of dust or dirt.
  - 1 Unplug the analyzer and remove the power cord from the back.
  - 2 Remove the air filter cover by pulling gently.
  - 3 Wash the filter in warm soapy water, rinse, and dry completely before reassembling.
  - 4 Replace clean and dry air filter into its cover and push the clips back into position.
  - 5 Plug the power cord back into the analyzer.
  - 6 Plug the power cord into the power source and press the Power button on the front. The analyzer will perform a diagnostic self test, which may take up to 4 minutes for the heaters to bring the disc chamber to operating temperature.

# Troubleshooting

The analyzer performs a series of internal quality checks to ensure accurate results. When a problem is detected, it either suppresses certain chemistry results (in which case it prints troubleshooting “flags” – such as ~~~, HEM, LIP, or ICT – in place of values) or cancels the run (the disc cancels and no results are printed). When either of these occurs, an error report can be printed or transmitted manually using the menu options.

The analyzer will display warning and error messages when problems occur. These messages include an internal error code that will assist Abaxis Technical Support in diagnosing the problem. Record the error message and/or print an error report before calling Abaxis Technical Support. The following are types of error codes and troubleshooting actions.

## Result Messages

- The symbols ~~~ are printed in place of numbers when a result cannot be determined. This may happen due to improper mixing of a reagent bead with diluted sample, a nonlinear reaction, an endpoint of a particular reaction not reached, a concentration outside the analyzer’s capabilities, or thermal stress. When a chemistry is suppressed (~~~), the analyzer will print an error report.
- **HEM, LIP, or ICT** is printed in place of the analyte concentration if the interferent has adversely affected the results. **LIP** is also printed if both lipemia and icterus have affected a particular analyte. **HEM** is also printed if hemolysis and icterus, hemolysis and lipemia, or hemolysis, lipemia, and icterus have affected a particular analyte. Examine the sample indices to determine if more than one interferent is affecting a particular result.
- **CONFIRM LOW RECOVERIES** with an exclamation mark (!) next to every analyte indicates that at least one analyte has a lower concentration than normally expected. If this occurs, re-run the sample. If it happens again, the sample may be problematic.

# Troubleshooting

## Electrostatic Discharge

An electrostatic discharge may cause the analyzer to “freeze up” while running a sample. If this happens, cancel the analysis manually, then power the analyzer off and back on again to restore operating conditions.

## Disc Cancellations

If the disc cancels, record the following information or print an error report and contact Abaxis Technical Support:

- Lot number
- Product name
- 4-digit error code
- Specimen type and sample

## Symbols on the Results Printout

The following table explains the symbols that may appear on the results printout. Call Abaxis Technical Support if the entire numerical line is displayed on the error report, including the line HEM (H), LIP (L), and ICT (I).

| Symbols printed next to the result |                                                                                                                                        |                                                                                                                                                                                                                                                                                                                           |
|------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Symbol                             | Meaning                                                                                                                                | Notes                                                                                                                                                                                                                                                                                                                     |
| *                                  | The result is outside the reference range.                                                                                             | An asterisk (*) is printed next to the analyte concentration if the results are outside the range. An asterisk is also printed if the result is at the upper or lower limit of the range, but because results are rounded, these are considered out-of-range as well. For example:<br>Na+ 145* (128–145) Cl- 98* (98–108) |
| <                                  | The result is lower than the dynamic range.                                                                                            | For more information on results with < or > printed before the value, print an error report. This report includes patient results, but the values approximate the analyte concentration.                                                                                                                                  |
| >                                  | The result is higher than the dynamic range.                                                                                           |                                                                                                                                                                                                                                                                                                                           |
| H                                  | After an LD value, indicates that hemolysis might affect the results.                                                                  |                                                                                                                                                                                                                                                                                                                           |
| < and H                            | If < appears before and H after an LD result, hemolysis is more extensive (100–150 mg/dL) and the true LD value is less than reported. |                                                                                                                                                                                                                                                                                                                           |
| C                                  | The result was calculated.                                                                                                             |                                                                                                                                                                                                                                                                                                                           |

### Symbols printed next to the result

|     |                                                                                                               |                                                                                                                         |
|-----|---------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------|
| !   | Confirm Low Recoveries. This indicates at least one analyte has a lower concentration than normally expected. | Re-run the sample. If the message recurs, results may fit patient profile clinically, or the sample may be problematic. |
| N/C | Not Calculated.                                                                                               |                                                                                                                         |

### Symbols (flags) printed instead of results

| Symbol     | Meaning                | Notes                                                                                                                                                                                                                                                                                                                                                                         |
|------------|------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| ~~~        | Chemistry interference | <b>HEM</b> , <b>LIP</b> , or <b>ICT</b> is printed in place of the analyte concentration if the interferent has adversely affected the results. <b>LIP</b> is also printed if both lipemia and icterus have been affected. <b>HEM</b> is also printed if hemolysis and icterus, hemolysis and lipemia, or hemolysis, lipemia, and icterus have affected a particular analyte. |
| <b>HEM</b> | Hemolysis interference |                                                                                                                                                                                                                                                                                                                                                                               |
| <b>LIP</b> | Lipemia interference   |                                                                                                                                                                                                                                                                                                                                                                               |
| <b>ICT</b> | Icterus interference   |                                                                                                                                                                                                                                                                                                                                                                               |

---

**ABAXIS USA**

3240 Whipple Rd.

Union City, CA 94587

Phone: 800-822-2947, option 2

Fax: 877-349-2087

Email: [ab-piccolotechsupport@zoetis.com](mailto:ab-piccolotechsupport@zoetis.com)

**ABAXIS Europe GmbH**

Bunsenstr. 9 – 11

64347 Griesheim – Germany

Phone: +49 6155 780 21 0

Support: +49 6155 780 21 300

E-Mail: [Int.PiccoloSupport@zoetis.com](mailto:Int.PiccoloSupport@zoetis.com)

## Support

---

Call ABAXIS Technical Support for questions regarding technical assistance or troubleshooting:

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**Abaxis Inc.**