



piccolo^{xpress}[®] chemistry analyzer

FAST, EASY, ACCURATE ON-SITE CHEMISTRY TESTING

REAL-TIME CHEMISTRY RESULTS FOR ONCOLOGY CENTERS

MINIMIZE THE TIME YOUR PATIENTS SPEND AT YOUR FACILITY

- Run chemistry results before and after treatment
- Up-to-the-minute information for titration
- Patients arrive later and go home earlier
- Reduce the risk of adverse events from chemo treatment
- Greater patient well-being and satisfaction

MAXIMIZE THROUGHPUT AND EFFICIENCY

- Minimize lost or erred results from the lab
- Ensure patient chemistries are always on time
- Reduce expensive chemotherapy drug waste
- Potential to see more patients per day
- Reduce staff burden of managing labs

○ FAST AND COMPLETE

From draw to complete panels in 3 easy steps and approximately 12 minutes



1 / Add sample



2 / Insert disc



3 / Read results

○ SIMPLE AND SMART

Fully automated system requires no special operating skills

- Self contained reagent discs deliver complete panel results in real time
- Very small sample – only 100 microliters of whole blood, serum or plasma
- As easy to use as a CD player – load sample, insert disc and read results
- Barcoded discs self calibrate with every run and minimize data input errors
- Built-in Intelligent Quality Control eliminates the need for daily external QC
- Results printed on sticky-backed roll-tape for easy placement in patient file
- 5 bi-directional USB ports for easy connectivity to EMR/LIS or peripherals
- Real-time analysis of electrolytes, kidney function, magnesium, LDH and more in a single panel

○ ACCURATE AND DEPENDABLE

Amazing technology ensures accurate results every time

- Precision and specificity comparable to expensive high-throughput analyzers
- Reference-grade enzymatic spectrophotometry
- Patented lyophilized reagent beads and space-age fluidics
- Absorption of reactions measured on 9 different wavelengths
- iQC[®] system provides continuous quality control checks on every run to eliminate the potential for erroneous results



3240 Whipple Rd., Union City, CA 94587 | 800.822.2947 | www.piccoloxpress.com

©2015 Abaxis, Inc., Piccolo Xpress[®] and iQC[®] are registered trademarks of Abaxis, Inc. 888-3024 Rev.D